

TRI-PMP-BXTM

Flow rate (per pump)	0.5ml to 105ml/min
Accuracy (per pump)	+/- 1%
Pump head (per pump)	8.1 meters (26.5')
Tubing size	Any 5mm O.D. tubing
Calibration	Single point
Modes of operation	Continuous dispensing Volume dispensing Dose over time Constant flow rate Dispense at startup
Data protocol	I²C
Default I ² C address	56 (0x38) 57 (0x39) 58 (0x3A)
Data format	ASCII
Operating voltage	12V 3 amp





Table of contents

EZO-PMP™ dimensions	3
Operating principle	4
Power consumption	4
Connectors	5
Connections to common CPUs	6
Sending commands to device	7
Requesting data from device	8
Response codes	9



I²C

Command quick page	10
Single report mode	11
Continuous dispensing	12
Volume dispensing	13
Dose over time	14
Constant flow rate	15
Dispense at startup	16
Pause dispensing	19
Stop dispensing	20
Invert dispensing direction	21
Total volume dispensed	22
Calibration	23
Enable/disable parameters	24
Device information	25
I ² C address change	26

Calibration theory	27
Accuracy	29
Viscosity	30
EZO-PMP™ tubing	31
Warranty	35



TRI-PMP-BX™ dimensions



I²C addresses


56 (0x38)

57 (0x39)

58 (0x3A)

Operating principle

The TRI-PMP-BX™ system consists of three individually controllable pumps that are connected together through a single I²C data line. It's best to think of this device as three individual smart pumps in one enclosure.



Power consumption

Logic	~40 mA (for all 3 pumps)
Motor	400 mA per pump when running


Absolute max ratings

Parameter	MIN	TYP	MAX
Storage temp (EZO-PMP™)	-65 °C		125 °C
Operational temp (EZO-PMP™)	-40 °C	25 °C	85 °C
Tubing life span	+1,000 hrs.		
Cassette life span	+1,500 hrs.		
Motor life span	+5,000 hrs.		

Connectors

The TRI-PMP-BX™ system requires two power sources:

- 3.3v – 5v logic voltage
- 12 – 24V motor voltage



I2C addresses **56 (0x38)** **56 (0x38)** **56 (0x38)**


Vcc **3.3V – 5.5V**

Clock speed **100 – 400 kHz**

SDA

SCL






0V **VCC** 0V



Data format

Readings	volume	Data type	floating point
Units	ml	Decimal places	2
Encoding	ASCII	Smallest string	3 characters
Format	string	Largest string	39 characters

Connections to common CPUs




3.3v	N/C
SDA	N/C
SCL	GND


Sending commands to device





Example



Advanced





Requesting data from device




413

SDA ← SDA ✓

SCL ← SCL ✓



Advanced




1 52 49 51 0 = 413

Dec ASCII Dec

Response codes

After a command has been issued, a 1 byte response code can be read in order to confirm that the command was processed successfully.

Reading back the response code is completely optional, and is not required for normal operation.



Example

```
I2C_start;  
I2C_address;  
I2C_write(EZO_command);  
I2C_stop;
```

```
delay(300);
```



```
Processing delay
```

```
I2C_start;  
I2C_address;  
Char[ ] = I2C_read;  
I2C_stop;
```

If there is no processing delay or the processing delay is too short, the response code will always be 254.

Response codes

Single byte, not string

255	no data to send
254	still processing, not ready
2	syntax error
1	successful request

Command quick reference

All commands are ASCII strings or single ASCII characters.

Command	Function	
Cal	performs calibration	pg. 23
D	dispense modes	pg. 12 – 18
i	device information	pg. 25
Invert	invert dispensing direction	pg. 21
I2C	change I ² C address	pg. 26
O	enable/disable parameters	pg. 24
P	pauses the pump during dispensing	pg. 19
R	returns a single reading	pg. 11
Tv	total volume dispensed	pg. 22
X	stop dispensing	pg. 20

Single report mode

Command syntax

300ms  processing delay

R returns a single value showing dispensed volume

Example

Response

R



1

Dec

2.50

ASCII

0

Null

(If issued half way through dispensing 5ml)



1

Dec

5.00

ASCII

0

Null

(If issued once dispensing has stopped)

Continuous dispensing

Pump on/pump off

300ms  processing delay

Command syntax

After running in continuous mode for 20 days the EZO-PMP™ will reset.

D,* dispense until the stop command is given

D,-* dispense in reverse until the stop command is given

D,? dispense status

Example

Response

D,*



1
Dec

0
Null

pump will continuously run at ~105ml/min (with supplied tubing)

D,-*



1
Dec

0
Null

pump will continuously run in reverse at ~105ml/min (with supplied tubing)

D,?



1
Dec

?D,10.00,1
ASCII

0
Null

Response breakdown

?D,* ,1

↑ ↑
last volume pump on
requested

Volume dispensing

Pump a specific volume

300ms  processing delay

Command syntax

where [ml] is any volume in millimeters ≥ 0.5

- D,[ml] dispense [this specific volume]
- D,[-ml] dispense [*in reverse* this specific volume]
- D,? dispense status

Example

Response

D,15



1 **0**
Dec Null

15 ml will be dispensed

D,-40.5



1 **0**
Dec Null

40.5 ml will be dispensed
in reverse

D,?



1 **?D,-40.50,0** **0**
Dec ASCII Null

Response breakdown

?D,-40.50,0

↑ last volume dispensed
↑ pump off

Dose over time

Pump a fixed volume over a fixed time

Command syntax

300ms  processing delay

D,[ml],[min] Dispense [this volume], [over this many minutes]

Example

Response


D,85,10


Wait 300ms

1
Dec

0
Null

Dispense 85ml over 10 mins



Constant flow rate

Maintain a constant flow rate

300ms  processing delay

Command syntax

After running in continuous mode for 20 days the EZO-PMP™ will reset.

DC,[ml/min], [min or *] [maintain this rate], [for this much time]

DC,? reports maximum possible flow rate

[ml/min] = a single number (int or float) representing the desired flow rate

[min or *] = the number of minutes to run or (*) indefinitely

A negative value for ml/min = reverse

Example

DC,25,40



1 0
Dec Null

Dispense 25ml per minute
for 40 minutes

DC,?




1 ?maxrate,58.5 0
Dec ASCII Null

The maximum flow rate is determined after calibration.
If the flowrate entered is too fast the EZO-PMP™ will send an error.

*TOOFAST

*ER



Dispense at startup

Pump a specific volume at startup and then stop

Use this command to make a simple fixed-volume pump

Command syntax

300ms  processing delay

Dstart,[ml] dispense [this specific volume] at startup

Dstart,off disables dispense at startup mode

Dstart,? startup dispense status

Example

Response


Dstart,10



Dstart,off



D,?



Continuous dispensing at startup

Pump on & continuously dispense

300ms  processing delay

Command syntax

After running in continuous mode for 20 days the EZO-PMP™ will reset.

- Dstart,*** dispense at startup until the stop command is given
- Dstart,-*** dispense in reverse at startup until the stop command is given
- Dstart,?** startup dispense status

Example

Response

Dstart,*



1 **0**
Dec Null

Pump will startup and continuously run at ~105ml/min (with supplied tubing)

Dstart,-*




1 **0**
Dec Null

Pump will startup and continuously run in reverse at ~105ml/min (with supplied tubing)

Dstart,?



1 **?Dstart,*** **0**
Dec ASCII Null



Dose Over time at startup

Pump a fixed volume over a fixed time at startup

Command syntax

300ms  processing delay

D,[ml],[min] Dispense [this volume], [over this many minutes] at startup

Example

Response

Dstart,85,10


Wait 300ms


1 **0**
Dec Null

Pump will startup and dispense 85ml over 10 minutes

Dstart,?


Wait 300ms

1 **?Dstart,85.00,10.00** **0**
Dec ASCII Null



Pause dispensing

300ms  processing delay

Command syntax

Issue the command again to resume dispensing

P pauses the pump during dispensing

P,? pause status

Example

P

 Wait 300ms
1 Dec 0 Null

P,?


 Wait 300ms
1 Dec ?P,1 0 Null
ASCII paused
or
 Wait 300ms
1 Dec ?P,0 0 Null
ASCII unpaused



dispensing



paused



dispensing

Stop dispensing

Command syntax

300ms  processing delay

X stop dispensing

Example

Response

X


Wait 300ms

1
Dec


***DONE,v**
ASCII

0
Null

v = volume dispensed



X
dispensing



***DONE,10.15**
Stopped

Invert dispensing direction

300ms  processing delay

Command syntax

Invert direction will be retained if power is cut

Invert changes dispensing direction of pump

Example


Response

Invert

 **Wait 300ms** **1** **0**
Dec Null

Invert,?

 **Wait 300ms** **1** **?Invert,1** **0** or  **Wait 300ms** **1** **?Invert,0** **0**
Dec ASCII Null Dec ASCII Null
inverted uninvverted



Total volume dispensed

Command syntax

300ms  processing delay

TV,? shows total volume dispensed

ATV,? absolute value of the total volume dispensed

Clear clears the total dispensed volume

Example

Response

TV,?

 **1** **?TV,623.00** **0**
Wait 300ms Dec ASCII Null

ATV,?

 **1** **?ATV,434.50** **0**
Wait 300ms Dec ASCII Null

clear

 **1** **0** **total now 0.00**
Wait 300ms Dec Null

This data will be lost if the power is cut.

Calibration

300ms  processing delay

Command syntax

Calibrate to the actual volume dispensed.

Cal,v v = corrected volume
Cal,clear delete calibration data
Cal,? device calibrated?

Example

Response



Cal,24.01




Wait 300ms **1** **0**
Dec Null

Cal,clear


Wait 300ms **1** **0**
Dec Null

Cal,?


Wait 300ms **1** **?Cal,1** **0** or 
Dec ASCII Null Dec ASCII Null
fixed volume volume/time


Wait 300ms **1** **?Cal,3** **0** or 
Dec ASCII Null Dec ASCII Null
both uncalibrated

Enable/disable parameters from output string

Command syntax


300ms  processing delay

O, [parameter],[1,0] enable or disable output parameter
O,? enabled parameter?

Example

Response

O,V,1

 **Wait 300ms** **1** **0**
Dec Null

enable volume being pumped

O,TV,0

 **Wait 300ms** **1** **0**
Dec Null

disable total volume pumped

O,ATV,1

 **Wait 300ms** **1** **0**
Dec Null

enable absolute volume pumped

O,?

 **Wait 300ms** **1** **?O,V,TV,ATV** **0**
Dec ASCII Null

if all three are enabled

Device information

Command syntax

300ms  processing delay

i device information

Example

Response

i



Wait 300ms

1

Dec

?i,PMP, 1.1

ASCII

0

Null

Response breakdown

?i, PMP, 1.1
↑ ↑
Device Firmware

I²C address change

Command syntax

300ms  **processing delay**

I2C,n sets I²C address and reboots into I²C mode

Example

Response

I2C,101

device reboot
(no response given)

Warning!

Changing the I²C address will prevent communication between the circuit and the CPU until the CPU is updated with the new I²C address.

Default I²C address are:

56 (0x38), **57 (0x39)**, **58 (0x3A)**


n = any number 1 – 127

Calibration theory

Uncalibrated accuracy +/- 5%

Calibrated accuracy +/- 1%

Before calibration is attempted all the air bubbles should be removed from the tubing. This is done by running the pump while tapping the tubing. If air bubbles are not removed from the tubing they will slowly group together into larger air bubbles. Over time this will lead to accuracy issues.




Calibration types

Volume calibration

Volume over time calibration


Calibration is optional. Both types of calibration are independent of each other and can be done at any time. Calibration can be done at any volume however; Atlas Scientific recommends using volumes above 5ml.

Equipment needed for calibration



An accurate graduated cylinder of at least 10ml.

Or



An accurate scale with a resolution of at least 0.1 grams

1 gram of water = 1ml
23.56 grams of water = 23.56ml

Calibration procedure


Calibration should be done with water and not a chemical

Make sure the tubing is full of water and has no bubbles before calibrating.

1. Instruct the pump to dispense a volume of water.
2. Measure the dispensed amount to determine how much water was actually dispensed.
3. Calibrate the pump by sending it the volume of liquid you have measured.

Example

Calibrate the pump by dispensing 10ml



1. Instruct the pump to dispense 10ml into a graduated cylinder or beaker on a scale.
2. Measure the amount of liquid that was actually dispensed.
3. Inform the pump how much liquid was actually dispensed.
4. Calibration is now complete.


Once the pump has been calibrated, it will accurately dispense any volume of liquid. Use the same procedure to perform a volume over time calibration.

Accuracy

Uncalibrated accuracy +/- 5%
Calibrated accuracy +/- 1%


Volume dispensing mode

calibrated at 10ml



Dose over time mode

calibrated at 10ml over 90 seconds




Viscosity

The EZO-PMP™ is capable of pumping liquids within a viscosity range of **0.1 – 2,000 cP**.


- 0.6 = Acetone
- 1 = Water
- 10 = Kerosene
- 100 = Corn Syrup
- 200 = Maple Syrup
- 2,000 = Honey
- 10,000 = Hershey Chocolate Syrup

Dispense accuracy



Pump head

Pump head refers to the maximum vertical height a pump can dispense. The EZO-PMP™ has a pump head of 8.1 meters (26.5').



EZO-PMP™ tubing

Tan tubing

Saint-Gobain™ PharMed™ BPT tubing

Length: 15.24cm

Outer diameter: 5mm

Inner diameter: 3mm

This tubing is highly chemically resistant and has 30X more resistant to mechanical wear than silicone tubing.



Inline tubing connectors

HDPE

Length: 2.54cm

Outer diameter: 8mm

Inner diameter: 2.8mm



Blue tubing

Silicone

Length: 2x 30.48cm


Outer diameter: 5mm

Inner diameter: 3mm


Bend radius: 15mm

Temperature -67°C to 200°C

Max pressure: 69 kPa (10 PSI)




Removing cassette




Removing tube assembly

The inner workings of the cassette are fragile and must be dismantled by hand. Using tools can damage or break the cassette.




Installing new tube assembly

1



Measure 75mm of pump tubing, and mark both ends with a soft-tip pen or marker.

2




Apply silicone lubricating grease to the marked areas on both the tubing and cassette axle.

Do not operate this device without lubrication!

Atlas Scientific recommends using **Super Lube** silicone lubricating grease.




3



Push the tube holder along with the greased and marked pump tube into the cassette until it clicks.


4



Gently pull out the pump tube, and insert the rotor into the pump tube. Align pump tube and rotor with the cassette axle.

Cassette

Roller




Make sure the pump tube is held between the roller and cassette.

✓ **Correct**



X Incorrect



Once the tubing has been replaced, run the pump for 3–5 minutes to break in the new tubing. **Remember, this pump can be run dry and does not need to pump liquid for the 3–5 minute break in period.**

Warranty

Atlas Scientific™ Warranties the TRI-PMP-BX™ to be free of defect during the debugging phase of device implementation, or 30 days after receiving the TRI-PMP-BX™ (which ever comes first).

The debugging phase

The debugging phase as defined by Atlas Scientific™ is the time period when the TRI-PMP-BX™ is inserted into a bread board, or shield. If the TRI-PMP-BX™ is being debugged on a bread board, the bread board must be devoid of other components. If the TRI-PMP-BX™ is being connected to a microcontroller, the microcontroller must be running code that has been designed to drive the TRI-PMP-BX™ exclusively and output the TRI-PMP-BX™ data as a serial string.

It is important for the embedded systems engineer to keep in mind that the following activities will void the TRI-PMP-BX™ warranty:

- Soldering any part of the TRI-PMP-BX™.
- Running any code, that does not exclusively drive the TRI-PMP-BX™ and output its data in a serial string.
- Embedding the TRI-PMP-BX™ into a custom made device.
- Removing any potting compound.

Reasoning behind this warranty

Because Atlas Scientific™ does not sell consumer electronics; once the device has been embedded into a custom made system, Atlas Scientific™ cannot possibly warranty the TRI-PMP-BX™, against the thousands of possible variables that may cause the TRI-PMP-BX™ to no longer function properly.

Please keep this in mind:

- 1. All Atlas Scientific™ devices have been designed to be embedded into a custom made system by you, the embedded systems engineer.**
- 2. All Atlas Scientific™ devices have been designed to run indefinitely without failure in the field.**
- 3. All Atlas Scientific™ devices can be soldered into place, however you do so at your own risk.**

Atlas Scientific™ is simply stating that once the device is being used in your application, Atlas Scientific can no longer take responsibility for the TRI-PMP-BX™ continued operation. This is because that would be equivalent to Atlas Scientific™ taking responsibility over the correct operation of your entire device.